



SK9000

*Double Chamber
Full-auto Hematology
Analyzer*



Features

- 3-part differentiation of WBC, 23 parameters.
- 8.4" color TFT screen, all testing parameter displayed simultaneously.
- Double channel counter, up to 60 samples test per hour.
- Sample volume : venous 9.8 ul, prediluted 20 ul for twice testing.
- Windows operation system, mouse and keyboard operation.
- Automatically sample probe cleannning (inside and outside).
- Large storage capacity: up to 10,000 samples +3 histograms.
- Internal thermal-sensitive printer or external printer.
- Electrical resistance for counting and SFT method of Hemoglobin.
- RS232 interface, PC connecting.





Main Technical Specification

- **Methodology:** Electrical resistance for counting, hemiglobincyanide method and SFT method for hemoglobin
- **Parameter:** 3-part differentiation of WBC; 20 parameters and 3 color histograms (WBC, RBC, PLT)
- **Work mode:** double Channel + unique Hemoglobin test system
Sample volume: 9.8 u L for Venous and capillary mode, 20 u L prediluted mode
- **Throughput:** More than 60 samples per hour, operatable 24 hours a day, auto sleeping and waking-up functions
- **Storage:** up to 100000 sample results including histograms can be stored, convenient for inquiry and management of history data
- **Operation language:** English
- **QC control:** X-B, L-J, X, SD, CV%
- **Reference Value Setting:** Male, Female, Children, Neonate
- **Input / Output:** RS232, parallel printer and keyboard
- **Print:** Graphic thermal printer with various printing format, optional external printer
- **Temperature:** 18 c – 30 c, wet <= 10-90 %
- **Power Supply:** 220 V + 22 VAC, 50+1HZ
- **Dimension:** 33 cm (L) * 38 cm (W) * 43 cm (H)
- **Weight:** 20 KG



Precision

Parameters	Range	Parameters	Range
WBC	0.0 -99.9x10 ⁹ /L	GRAN#	0 - 99.9x10 ⁹ /L
RBC	0.00- 9.99x10 ¹² /L	HCT	0.0 - 100.0 %
HGB	0.00 - 300g/L	MCH	0.0 - 999.9pg
PLT	0 - 3000x10 ⁹ /L	MCHC	0.0 - 999.9g/L
MCV	0 - 250fL	RDW-SD	0.0 - 99.9 fL
LYM%	0 - 100 %	RDW-CV	0.0 - 99.9 %
MID%	0 - 100 %	PDW	0.0 - 30.0 %
GRAN%	0 - 100 %	MPV	0.0 - 30.0fL
LYM#	0 - 99.9x10 ⁹ /L	PCT	0.0 - 9.99 %
MID#	0 - 99.9x10 ⁹ /L	P-LCR	0.0 - 99.9 %



Parameters

WBC, Lymph#, Mid#, Gran#, Lymph%, Mid%, Gran%, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-CV, RDW-SD, PLT, MPV, PDW, PCT, L-PCR
3- histograms : WBC, RBC and PLT